Liquid nitrogen storage system

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
Liquid Nitrogen Biological Container Monitor System

- The control system adopts the control mode of customized PLC + touch screen, which is easy to operate and the system is stable and reliable.
- Two control modes can be carried out by the system: manual control, automatic control.
- High-precision monitor based on differential pressure system Unique design to estimate and track liquid nitrogen consumption.
- Up to 17 audio/visual alarm systems.
- The hot gas path maintains a stable temperature gradient and reduces excess liquid nitrogen evaporation.
- Can store 30,000 usage records with ten years of storage capacity.

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STAINLESS STEEL TANK SERIES

The new wide-neck stainless steel liquid nitrogen bio-tank provides excellent protection for biological samples. The microprocessor-based temperature and liquid nitrogen level monitor system closely monitors equipment operation to ensure sample safety.

The standard integrated folding steps, auxiliary workbench and ergonomic cover design ensure the safety and simplicity of the equipment. Compatible with liquid nitrogen gas and liquid phase storage mode, ensuring excellent temperature uniformity and stability throughout the entire cabinet even under extended opening conditions.

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Humanization design, comfortable operation lockable cover, fully guarantee the safety of the samples

- The hot gas bypass function automatically discharges nitrogen from the pipeline when liquid nitrogen is injected, which reduces the influence on the internal temperature of the container and ensures the safety of the stored samples.
- A wide range of capacities to choose from, maximizing your storage needs. Store more samples per unit area, truly achieve high-density storage and reduce storage costs.
- Excellent temperature uniformity and stability, gas storage top temperature is less than -196 degrees Celsius. Excellent temperature stability after opening locking cover, top temperature is less than -160 degrees Celsius after opening lacking cover for 4 hours.
- Advanced temperature, liquid level monitor and alarm system, and network monitor to automatically replenish liquid nitrogen, safe and
- labor-saving One-button defogging function makes it easy to find and pick up samples.

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The Cryobio 400 series is a small-scale high-density storage, providing users with fully automatic, safe and reliable ultra-low temperature liquid nitrogen storage system, which can store 15600~19500 2ml cryo-tubes. When the large volume of liquid nitrogen at the bottom of the tray is designed for gas phase storage, the temperature difference of the entire storage area does not exceed 10 degrees Celsius, and the lowest temperature of the top of the racks can reach -190 degrees Celsius.

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PRODUCT FEATURES

- Automatic rehydration system and hot gas discharge system.
- Standard vacuum inlet hose.
- Gas phase storage, top up to -190°C.
- Optional removable single or double step.
- Optional battery backup.
- Five-year vacuum warranty.

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Static evaporation and static hold time are theoretical values, actual evaporation and retention time will be affected by usage, atmospheric conditions and manufacturing tolerances

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The Cryobio 800 series is a medium-scale high-density storage, providing users with fully automatic, safe and reliable ultra-low temperature liquid nitrogen storage system, which can store 36400~42000 2ml cryo-tubes. When the large volume of liquid nitrogen at the bottom of the tray is designed for gas phase storage, the temperature difference of the entire storage area does not exceed 10 degrees Celsius, and the lowest temperature of the top of the racks can reach -190 degrees Celsius.

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PRODUCT FEATURES

- Automatic rehydration system and hot gas discharge system.
- Standard vacuum inlet hose.
- Gas phase storage, top up to -190℃.
- Standard double step.
- Optional battery backup.
- Five-year vacuum warranty.

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
The Cryobio 1800 series is a medium-scale high-density storage, providing users with fully automatic, safe and reliable ultra-low temperature liquid nitrogen storage system, which can store 79950~94500 2ml cryo-tubes. When the large volume of liquid nitrogen at the bottom of the tray is designed for gas phase storage, the temperature difference of the entire storage area does not exceed 10 degrees Celsius, and the lowest temperature of the top of the racks can reach -190 degrees Celsius.

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PRODUCT FEATURES

- Automatic rehydration system and hot gas discharge system.
- Standard vacuum inlet hose.
- Gas phase storage, top up to -190°C.
- Standard double step.
- Optional battery backup.
- Five-year vacuum warranty.

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Aluminum alloy container is composed of outer shell, inner vessel, neck tube, multi-layer thermal isolation, racks etc. The outer shell of the container is made of aluminum alloy plate, which is light weight, high strength at low temperature and corrosion resistance. The neck tube is made of fiberglass which is high mechanical strength and small thermal conductivity.

The multi-layer insulator adopts an aluminum foil with particularly excellent reflection performance as a reflective screen, and a material having a small thermal conductivity and a low gas release rate is used as a heat insulating material to reduce heat radiation.

The interlayer between the outer shell and the inner tank is in a high vacuum state to prevent the heat convection of the gas, and the adsorbent having a large adsorption amount at a low temperature is used to ensure the long-term stability and reliability of the performance of the container, and the product life is not less than five years.

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WIDE-NECK RACK SERIES

Wide-neck rack series combines the advantages of low liquid nitrogen consumption with medium-range capacity cryo-tube. Larger neck diameter and easy-to-extract racks are convenient for long-term storage of biological samples and frequent sample extraction.

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PRODUCT FEATURES

- Wider neck Large capacity low liquid nitrogen consumption.
- Safety locking cover.
- High strength, lightweight aluminum structure.
- Five-year vacuum warranty.

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THE STATIC LARGE-CAPACITY STORAGE SERIES

The static large-capacity series is specially designed for samples that require long-term static storage. There are two types of products with large capacity or long storage period. The products are durable, light weight, world-leading manufacturing technology, and have five years of vacuum warranty and many other advantages.

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
PRODUCT FEATURES

- Designed for large capacity.
- Ultra low evaporation loss.
- Convenient and light packing.
- Large capacity low liquid nitrogen consumption Safety locking cover.
- High strength, lightweight aluminum structure.
- Five-year vacuum warranty.

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The Static Large-capacity Storage Series													
Models	YDS-10-80	YDS-10-125	YDS-13	YDS-15	YDS-15-125	YDS-20	YDS-25	YDS-30	YDS-30-80	YDS-35	YDS-35-80	YDS-50	YDS-50-80
Maximum Storage Capacity													
Numbers of Containers	6	-	6	6	-	6	6	6	6	6	6	6	6
Storage Working Capacity (0.5ml)	2244	-	1284	792	-	1284	1284	1284	3624	1284	3624	1284	3624
Storage Working Capacity (0.25ml)	5022	-	2832	1788	-	2832	2832	2832	8460	2832	8460	2832	8460
Performance													
Volume L	10	10	13	15	15	20	25	30	30	35	35	50	50
Static evaporation rate (L/day)	0.21	0.42	0.12	0.11	0.36	0.12	0.12	0.12	0.21	0.12	0.22	0.23	0.34
Static Liquid Nitrogen Safe Period (days)	48	24	109	134	42	168	208	254	147	286	159	213	147
Dimensions													
Inlet Diameter (mm)	80	125	50	50	125	50	50	50	80	50	80	50	80
Total Height (mm)	578	578	618	594	594	675	700	648	651	672	675	744	747
Outer Diameter (mm)	300	300	300	396	396	396	396	451	451	451	451	508	508
Container Height (mm)	120	-	276	120	-	276	376	276	276	276	276	276	276
Container Diameter (mm)	38	-	38	38	-	38	38	38	63	38	63	38	63
Weight (kg)	5.9	5.9	6.3	8.5	-	9.6	10.7	13	13.2	14.2	14.5	15.4	15.7

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THE PORTABLE SMALL-CAPACITY SERIES


The portable small-capacity series is designed for small capacity users, mainly used for portable transportation of livestock, frozen sperm breeding storage and biological samples.

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PRODUCT FEATURES

- Suitable for small capacity liquid storage.
- Low liquid nitrogen evaporability and easy maintenance.
- Standard safety locking cover.
- High strength, low weight aluminum construction.
- Five-year vacuum warranty.

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The Portable Small-capacity Series						
Models	YDS 2 30	YDS 2 35	YDS 3	YDS 6	YDS 8 80	YDS 10
Maximum Storage Capacity						
Numbers of Canisters	3	3	6	6	6	6
Straw Holding Capacity 0.5ml	90	165	792	792	2244	792
Straw Holding Capacity 0.25ml	204	330	1788	1788	5022	1788
Performance						
Volume L	2	2	3	6	8	10
Static evaporation rate (L/D)	0.07	0.08	0.12	0.12	0.21	0.12
Static Liquid Nitrogen Hold Time(D)	28	24	26	52	38	86
Dimensions						
Neck Diameter (mm)	30	35	50	50	80	50
Total Height (mm)	402	402	427	468	498	528
Outer Diameter (mm)	197	197	223	300	300	300
Canister Height (mm)	120	120	120	120	120	120
Canister Diameter (mm)	19	25	38	38	63	38
Empty Weight (mm)	3.2	3.2	4.4	5.4	7.2	6.6

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THE AIR TRANSPORTATION SERIES

The air transportation series is suitable for transporting samples in cryogenic environment (gas phase storage, temperature below -190°C), which avoids the danger of liquid nitrogen release during transportation, especially for short-time air transportation such as airplanes. There is liquid nitrogen adsorbent inside, which can absorb and retain liquid nitrogen.

Even if the container is dumped, there will be no liquid nitrogen overflow.

There is also a special stainless steel mesh screen to separate the storage space and the liquid nitrogen absorbent, which can avoid the liquid nitrogen absorbent material mixed into the cryo-tubes. It is mainly suitable for user needs of laboratory users and a small number of samples for short-term transportation.

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PRODUCT FEATURES

- Gas phase cryogenic storage.
- Fast liquid nitrogen filling.
- High strength aluminum structure.
- Lockable cover.
- Five-year vacuum warranty.

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LIQUID NITROGEN STORAGE AND DISTRIBUTION SERIES

The liquid nitrogen storage and distribution series is designed for small liquid nitrogen storage and short-distance transportation. It has small neck and low evaporation loss rate which is an efficient, convenient and economical way to store and distribute liquid nitrogen. It can be dumped during operation. Pressing liquid nitrogen. Liquid nitrogen can be discharged by pouring and pressurizing during operation.

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PRODUCT FEATURES

- Low liquid nitrogen consumption.
- Convenient and light weight packaging.
- Large capacity low liquid nitrogen consumption.
- High strength, lightweight aluminum construction.
- Five-year vacuum warranty.

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